 CRESPINE® GEL

CRESPINE® GEL

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BP_CrespineGel_BPZ_003

CRESPINE® GEL

Injectable implant for the treatment of osteoarthritis

Hyaluronic Acid Implant of High Purity – Single Injection

Instructions for use

Description

CRESPINE®GEL is an absorbable intra-articular implant of high purity. It is a medical device, manufactured from hyaluronic acid of non-animal origin.

CRESPINE®GEL is a sterile, apyrogenic, visco-elastic, water-insoluble, biologically compatible (non-immunising, non-inflammatory, non-toxic) gel implant, manufactured from cross-linked hyaluronic acid that is obtained by fermentation. Hyaluronic acid is a naturally occurring polysaccharide in the dermal matrix of the human skin and a component of synovial fluid. Hyaluronic acid is chemically, physically and biologically identical in the tissues of all higher organisms.

Composition

1 ml CRESPINE®GEL contains:

Sodium hyaluronate, crosslinked	14.0 mg
Sodium chlorid	6.9 mg
Water for injections ad	1.0 ml

Mode of action

CRESPINE®GEL is injected into joints as a visco-supplement for the reduction of pain symptoms and improved functioning of the joints.

Indication and application

CRESPINE®GEL is indicated for pain and joint function disorders as a result of degenerative deficits as a single injection treatment course. After the injection hyaluronic acid forms a layer on the surface of the joint cartilage and improves the pain symptoms.

Treatment areas

Osteoarthritis of knee or hip

Contraindications

CRESPINE®GEL must not be used in the event of:

- autoimmune illnesses
- incompatibility with gram-positive bacteria
- active inflammation or infection processes
- anticoagulant therapies
- known allergy to hyaluronic acid

Precautions for use

CRESPINE®GEL is only intended for intra-articular injection and must not be injected into blood vessels. This could lead to an occlusion of the vessels and embolism.

Clinical data concerning breastfeeding women and young persons under 18 years of age are not available. This fact should generally be taken into account before administering **CRESPINE®GEL** to women of childbearing potential, especially during early pregnancy when maximum organogenesis takes place.

CRESPINE®GEL is supplied in a sterile syringe for one single treatment and is ready for use. If the packaging is opened/defective, it must not be used.

It must not be mixed with other injectable agents.

Precautions

As with all injection treatments of the joints, the general, valid contraindications for intra-articular injection are also to be observed in treatment with **CRESPINE®GEL**.

The skin must be disinfected prior to treatment. As with all implants, **CRESPINE®GEL** must not be used in patients with existing infections and inflammation processes in the vicinity of the implantation site. Patients should not take any aspirin, steroids or high doses of vitamin E before treatment, as these substances can lead to bleeding and a tendency to inflammation at the injection site. The syringes and needles used are considered contaminated and must therefore be destroyed in accordance with the recognised rules of medical practice.

Effect on the ability to drive and use machines:

It is up to the physician to decide whether the patient is able to drive and to use machines, depending on the patient's general state of health after application.

Adverse effects

Caused by the injection:

As with any other injection, patients may suffer from the following symptoms:

- temporary erythema
- slight swelling
- pain
- sensation of heat

Typically these reactions spontaneously disappear within 2 to 5 days following the injection.

Caused by the medical device:

Hypersensitivity to hyaluronic acid after the injection is reported in less than 1 %. This hypersensitivity manifests itself by extensive erythema and swelling. These reactions can occur immediately after the injection or up to 2 to 4 weeks later. The clinical data show that these hypersensitivity reactions are mild or moderate and last for a maximum of two weeks. Patients with multiple allergies should be excluded from treatment. Hypersensitivity reactions can be caused by endotoxin residues

in the hyaluronic acid solution. Therefore the hyaluronic acid used in **CRESPINE®GEL** is endotoxin free (approved by independent laboratory, GLP-certified) and possesses the highest purity for hyaluronic acid products.

Application instructions

A full anamnesis must be performed before treatment, in order to rule out any contraindications. The practitioner has to advise the patient of all precautions and possible adverse effects prior to treatment. The area to be treated must be carefully prepared with an antiseptic before treatment. The syringe is removed from the blister pack, the cap covering the tip of the syringe is removed from the syringe and a needle is applied to the syringe.

CRESPINE®GEL is wrapped in sterile packaging and is intended for single use only. It must not be mixed with other injection agents. Do not re-sterilise **CRESPINE®GEL**. The needles used for implantation have to be selected according to their suitability for intra-articular injection.

To minimize the likelihood of intravascular injection, aspiration should be performed before the application. If blood is aspirated, the needle must be repositioned until no reflow of blood can be elicited by aspiration. Please note, however, that the absence of blood in the syringe does not ensure the avoidance of intravascular injection.

Advice

The injection technique is of decisive importance for the final outcome.

CRESPINE®GEL must only be used by authorised medical personnel in agreement with the regional legal regulations. It is possible to treat several joints at the same time. The common contraindications for intra-articular injection have to be observed.

CRESPINE®GEL must not be re-sterilised. If the packaging is open/defective or if the blister is damaged, it must not be used.

Important information

The graduation on the syringe is an orientation aid for the user and refers to the final volume. It does not possess a measuring function, but instead only states the quantities used, relating to the nominal volume of 2 ml.

Nature and contents of container

CRESPINE®GEL is supplied in sterile packaging in a 2 ml syringe with integrated Luer lock adapter in a blister, for single use only. Instructions for use and labels with the batch number and date of use are packed together with the blister in a cardboard box. The patient is given one of the labels to guarantee traceability of the product.

Precautions for storage

Do not freeze or expose to extreme heat.

CRESPINE®GEL should be stored at room temperature 18 °C – 25 °C.

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BioPolymer GmbH & Co. KG